

# Isolation Gown, SS



Product name	: Isolation Gown (nonsterile)
Fabric material	: SS Non-woven material
Fabric weight	: 20gsm, 25gsm, 30gsm, 35gsm, 40 gsm
Sleeve	: Short or Long
Back neck	: Velcro and/or tie upon request
Sleeves	: Knitted cuffed sleeves or elastic at the cuffs
Style	: Latex free elastic cuffs or knitted cuffs
Colors	: White, blue, green
Size	: S, M, L, XL, 2XL, 3XL. Custom made size available
Features	: Antistatic, Fluid repellend, metal free
Useage	: For medical or non medical use, food processing, nursing
Latex	: Latex Free, ideal for people with potential latex reactions
Certificate	: CE, ISO 9001:2015, ISO 13485:2016
Reference	: G5.06(20gsm), G5.07(25gsm), G5.08(30gsm), G5.09(35gsm), G5.10(40gsm)
Weight	: Depends on size and gsm of material
Package	: Carton 40x60x40 (LxWxH) 80-120 pcs/carton depends on size and gsm

Size	A	B	Tolerance	Ref.Code 20 gsm	Ref.Code 25 gsm	Ref.Code 30 gsm	Ref.Code 35 gsm	Ref.Code 40 gsm
S	110	135	± 2	G5.06.S	G5.07.S	G5.08.S	G5.09.S	G5.10.S
M	120	140	± 2	G5.06.M	G5.07.M	G5.08.M	G5.09.M	G5.10.M
L	125	145	± 2	G5.06.L	G5.07.L	G5.08.L	G5.09.L	G5.10.L
XL	130	150	± 2	G5.06.XL	G5.07.XL	G5.08.XL	G5.09.XL	G5.10.XL
2XL	135	155	± 2	G5.06.2XL	G5.07.2XL	G5.08.2XL	G5.09.2XL	G5.10.2XL
3XL	145	160	± 2	G5.06.3XL	G5.07.3XL	G5.08.3XL	G5.09.3XL	G5.10.3XL

