



Type 5 & 6

Product name	: Isolation Coverall, Laminated, nonsterile
Fabric material	: PP quality spunbond laminated with PE (PP+PE)
Fabric weight	: 40-45gsm or 55-57gsm
Weight	: Depends on size and gsm of material
Colors	: White
Size	: As per brlow table. Custom made size available
Features	: Antistatic, waterproof, impermeable to microbial, metal free
Useage	: For medical or non medical use, food processing, nursing
Protection	: Excellent barrier against dry particles and liquid chemical splash
Comfort	: Elasticated or knitted cuffs and elasticated waist and ankles for increased comfort and freedom of movement
Zipper	: Zipper with fully sealable storm flap
Strip	: Optional strip or stripless production
Latex	: Latex Free, ideal for people with potential latex reactions
Certificate	: ISO 9001:2015, ISO 13485:2016, CE, EN 14126, EN 14325
Package	: Carton 40x60x40 (LxWxH) 50-60 pcs/carton depends on size and gsm

Size	Height cm	Chest Girth cm	Tolerance	Ref.Code 40-45 gsm	Ref.Code 55-57 gsm
XS	156-164	76-84	± 2	C3.01.XS	C3.02.XS
S	162-170	84-92	± 2	C3.01.S	C3.02.S
M	168-176	92-100	± 2	C3.01.M	C3.02.M
L	174-182	100-108	± 2	C3.01.L	C3.02.L
XL	180-188	108-116	± 2	C3.01.XL	C3.02.XL
2XL	186-194	116-124	± 2	C3.01.2XL	C3.02.2XL
3XL	192-200	124-132	± 2	C3.01.3XL	C3.02.3XL
4XL	200-208	132-140	± 2	C3.01.4XL	C3.02.4XL



EN 13982 (Type 5B)
Protection against
solid particles



EN 13034 (Type 6B)
Protection against
liquid chemicals



EN 1149 -5
Anti Static



EN 14126
Protection against
infective agents

