

# Disposable Face Mask



Product Name	: Disposable Face Mask
Condition	: Non-sterile
Level	: Type I, Type II, Type IIR
BFE	: Type I ( $\geq 95\%$ ), Type II & IIR ( $\geq 98\%$ )
Colours (outer layer)	: White, light blue, light green or customer defined
Material, middle layer	: Melt-blown or Crownblond for filtering
Material, outer layer	: Non-woven spunbond
Material, inner layer	: Non-woven spunbond
Earloop	: 3mm latex-free round elastic earloop
Tie behind	: Ultrasonically bonded spunbond
Nose-piece	: Plastic covered metal
Fixing method	: Full ultrasonic, seamless
Dimensions (body)	: 170mm x 95mm (L x H)
Test Method	: EN 14683
Certificates	: CE, ISO 9001:2015, ISO 13485:2016

Reference Codes	Type I	Type II	Type IIR
Earloop	M0.01	M0.02	M0.03
Tie behind	M0.11	M0.12	M0.13