



Technical Specifications

Product name	: Disposable PE Glove
Glove material	: Polyethylene (HDPE)
Thickness	: 15-20 micron
Color	: Transparent
Features	: Easy-to-use, powder-free, latex free, economic product. It is used wherever hygiene is important.
Right-Left hand	: Compatible with both hands
Bag	: 100 pcs in a plastic pouch
Cartoon	: 10,000 pcs (100 plastic pouch)
Cartoon dimension	: 28x38x21.5 cm
Cartoon weight	: 7.5 kg gross weight (7.0 kg net weight)
Certificate	: CE, ISO 9001:2015, ISO 13485:2016