



Product name	: Disposable PE Apron, Without Sleeve
Material	: PE (LDPE)
Thickness	: 32 micron or 16 micron
Weight	: 25 gr \pm %5 (32 mic.), 12.5 gr \pm %5 (16 mic.),
Style	: Without sleeve
Colors	: Transparent, white, blue
Size	: L
Water isolation	: Waterproof 100%
Features	: Isolate water, paints, blood, alcohol, bacteria, virus
Useage	: For medical or non medical use, food processing, nursing
Latex	: Latex Free, ideal for people with potential latex reactions
AAMI Level	: Level 3 (Only on the apron)
Certificate	: CE, EN 811, ISO 9001:2015, ISO 13485:2016
Product Code	: A0.01 : 32 micron A0.02 : 16 micron