

Disposable PE Gown - Apron



Product name	: Disposable PE Gown, Apron
Material	: PE (LDPE)
Thickness	: 30 micron
Weight	: 55 gr
Sleeve	: Long
Style	: Thumbhole sleeves to avoid glove "pull-out"
Colors	: Transparent, white, blue
Size	: M, L, XL
Water isolation	: Waterproof 100%
Features	: Isolate water, paints, blood, alcohol, bacteria, virus
Useage	: For medical or non medical use, food processing, nursing
Latex	: Latex Free, ideal for people with potential latex reactions
AAMI Level	: Level 4
Certificate	: CE, EN 811, ISO 9001:2015, ISO 13485:2016
Product Code	: G4.01