




Product name	: Disposable Gown, Laminated	<b>Type 5 &amp; 6 Cat III</b>
Fabric material	: PP quality spunbond laminated with PE (PP+PE)	
Fabric weight	: 40-45gsm or 55-57gsm	
Weight	: Depends on size and gsm of material	
Sleeve	: Long	
Back neck	: Velcro or tie upon request	
Style	: Latex free elastic cuffs or knitted cuffs	
Colors	: White, blue	
Size	: S, M, L, XL, 2XL, 3XL. Custom made size available	
Features	: Antistatic, waterproof, impermeable to microbial, metal free	
Useage	: For medical or non medical use, food processing, nursing	
Latex	: Latex Free, ideal for people with potential latex reactions	
AAMI Level	: Level 4	
Certificate	: CE, EN 811, ISO 9001:2015, ISO 13485:2016	
Reference	: G3.04 (40-45gsm), G3.05 (55-57gsm)	
Package	: Carton 40x60x40 (LxWxH) 50-60 pcs/carton depends on size and gsm	

Size	A	B	Tolerance	Ref.Code 40-45 gsm	Ref.Code 55-57 gsm
S	110	135	± 2	G3.04.S	G3.05.S
M	120	140	± 2	G3.04.M	G3.05.M
L	125	145	± 2	G3.04.L	G3.05.L
XL	130	150	± 2	G3.04.XL	G3.05.XL
2XL	135	155	± 2	G3.04.2XL	G3.05.2XL
3XL	145	160	± 2	G3.04.3XL	G3.05.3XL



EN 13982 (Type 5B)  
Protection against  
solid particles



EN 13034 (Type 6B)  
Protection against  
liquid chemicals



EN 1149 -5  
Anti Static



EN 14126  
Protection against  
infective agents