



- Product name : Isolation Gown (nonsterile)
 Fabric material : SMS Non-woven material
 Fabric weight : 30gsm, 35gsm, 40gsm SMS
 Sleeve : Long
 Back neck : Velcro and/or tie upon request
 Sleeves : Knitted cuffed sleeves or elastic at the cuffs
 Style : Latex free elastic cuffs or knitted cuffs
 Colors : White, blue, green
 Size : S, M, L, XL, 2XL, 3XL. Custom made size available
 Stitch Method : Ultrasonically stitched or Thread stitched
 Features : Antistatic, Fluid repellend, metal free
 Useage : For medical or non medical use, food processing, nursing
 Latex : Latex Free, ideal for people with potential latex reactions
 Certificate : CE, ISO 9001:2015, ISO 13485:2016
 Reference : Thread Stitch : G0.02(30gsm), G0.03(35gsm), G0.04(40gsm)
 Ultrasonic Stitch : G0.08(30gsm), G0.09(35gsm), G0.10(40gsm)
 Weight : Depends on size and gsm of material
 Package : Carton 40x60x40 (LxWxH) 70-90 pcs/carton depends on size and gsm

Size	A	B	Tolerance	Thread Stitch			Ultrasonic Stitch		
				Ref.Code 30 gsm	Ref.Code 35 gsm	Ref.Code 40 gsm	Ref.Code 30 gsm	Ref.Code 35 gsm	Ref.Code 40 gsm
S	110	135	± 2	G0.02.S	G0.03.S	G0.04.S	G0.08.S	G0.09.S	G0.10.S
M	120	140	± 2	G0.02.M	G0.03.M	G0.04.M	G0.08.M	G0.09.M	G0.10.M
L	125	145	± 2	G0.02.L	G0.03.L	G0.04.L	G0.08.L	G0.09.L	G0.10.L
XL	130	150	± 2	G0.02.XL	G0.03.XL	G0.04.XL	G0.08.XL	G0.09.XL	G0.10.XL
2XL	135	155	± 2	G0.02.2XL	G0.03.2XL	G0.04.2XL	G0.08.2XL	G0.09.2XL	G0.10.2XL
3XL	145	160	± 2	G0.02.3XL	G0.03.3XL	G0.04.3XL	G0.08.3XL	G0.09.3XL	G0.10.3XL

