

# Particulate Respirator



Product Name	: Particulate Respirator
Condition	: Non-sterile, disposable
Type	: FFP1, FFP2, FFP3
Filtration Efficiency	: $\geq 80\%$ (0.3 microns) FFP1 $\geq 94\%$ (0.3 microns) FFP2 $\geq 95\%$ (0.3 microns) FFP3
Colours (outer layer)	: White, light blue, light green
Straps	: Featured latex free knitted
Headbands or Ear Loops	: Optional
Exhalation Valve	: Optional
Nose-piece	: Zinc and iron alloy or Aluminium
Fixing method	: Full ultrasonic, seamless
Test Method	: EN 149
Certificates	: CE, ISO 9001:2015, ISO 13485:2016