



Product names	: Shoe Cover, Nonwoven, standart Shoe Cover, Nonwoven, with non-slip sole
Material	: PP (SS)
Thickness	: Standart : 25 gsm SS. Non-slip sole : 25 gsm SS for sides, 40 gsm for sole Different fabrics can be applied upon request
Elastic	: Double elastic on the ankles
Sewing	: Standart : Ultrasonic Non-slip sole : Overlock sewing
Type of non-slip	: Non-slip rubber sole
Colors	: White, blue, green
Size	: 50 cm single size
Water isolation	: Waterrepellent
Useage	: For medical or non medical use, food processing, nursing
Latex	: Latex Free
Certificate	: ISO 9001:2015, ISO 13485:2016
Packing	: 800 pcs in a carton 40x40x40cm. (8 pouches, 100 pcs in a poche)
Carton Weight	: 3.20 kg for standart 4.80 kg for non-slip sole
Product Code	: SC.01 : Shoe Cover, Nonwoven, standart SC.02 : Shoe Cover, Nonwoven, with non-slip sole