


# Reinforced Surgical Gown



- Product name : Reinforced Surgical Gown (Sterile or Nonsterile)  
 Fabric material : SMS Non-woven material + PE (Polyethylene) Film Reinforced  
 Fabric weight : 40gsm SMS + 18-20gsm PE  
 Sleeve : Long  
 Back neck : Velcro and Tie  
 Sleeves : Knitted cuffed sleeves  
 Colors : White, blue, green  
 Size : S, M, L, XL, 2XL, 3XL. Custom made size available  
 Stitch Method : Ultrasonically stitched or Thread stitched  
 Features : Antistatic, Fluid repellend, metal free  
 Useage : For medical or non medical use, food processing, nursing  
 Latex : Latex Free, ideal for people with potential latex reactions  
 Certificate : CE, ISO 9001:2015, ISO 13485:2016  
 Reference : Thread Stitch : G6.04  
                   Ultrasonically Stitch : G6.05  
 Weight : Depends on size and gsm of material  
 Package : Carton 40x60x40 (LxWxH) 50-60 pcs/carton depends on size

Size	A	B	Tolerance	Thread Stitch	Ultrasonic Stitch
				40gsm SMS + PE	40gsm SMS + PE
				Ref. Code	Ref. Code
S	110	135	± 2	G6.04.S	G6.05.S
M	120	140	± 2	G6.04.M	G6.05.M
L	125	145	± 2	G6.04.L	G6.05.L
XL	130	150	± 2	G6.04.XL	G6.05.XL
2XL	135	155	± 2	G6.04.2XL	G6.05.2XL
3XL	145	160	± 2	G6.04.3XL	G6.05.3XL

